UrinAX CL/MC/E

Transportmedium, Keimzahlbestimmung und selektive Anzucht von Bakterien im Urin
Milieu de transport, dosage des germes et développement sélectif de bactéries dans l'urine
Mezzo di trasporto, conta batterica e moltiplicazione selettiva dei batteri nell'urina
Transport medium, bacterial count and selective propagation of bacteria in urine

Gebrauchsanweisung Notice d'utilisation Istruzioni per l'uso Instructions for use

REF AL0203
BD257499 Cont.

- 10 Dip-Slides / Lames immergées /
Piastrine a immersione / Dip slides
- 10 Patientenetiketten / Etiquettes pour patient /
Etichette paziente / Patient labels

In vitro diagnosticum
Usage in vitro
Per uso diagnostico in vitro
For in vitro diagnostic use

Axon Lab AG, CH-5405 Baden
Tel. +41 (56) 484 80 80
Fax +41 (56) 484 80 99
info@axonlab.ch · www.axonlab.ch

Schweiz
Axon Lab AG
Täferstrasse 15
CH-5405 Baden
Tel. +41 (0)56 484 80 80
Fax +41 (0)56 484 80 99
info@axonlab.ch
www.axonlab.ch

Suisse
Axon Lab SA
En Budron E9
CH-1052 Le Mont-sur-Lausanne
Tél. +41 (0)21 651 71 50
Fax +41 (0)21 651 71 59
info_f@axonlab.ch
www.axonlab.ch

Deutschland
Axon Lab AG
Heinrich-Otto-Strasse 1
D-73262 Reichenbach/Stuttgart
Tel. +49 (0)7153 9226-0
Fax +49 (0)7153 9226-12
info@axonlab.de
www.axonlab.de

Österreich
Axon Lab AG
Gewerbepark 1
A-6404 Polling
Tel. +43 (0)523 88 77 66
Fax +43 (0)523 88 77 66-12
info@axonlab.at
www.axonlab.at

UrinAX CL/MC/E

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Bibliografia

NCCLS Publication: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Volume 24:14, 2004.

Abkürzungen/Abréviations/Abbreviazioni/Abbreviations

Symbole/Symboles/Simboli/Symbols

Symbol/Symbole/Simbolo	Bedeutung/Signification/Significato
	Bestellnummer/Numéro d'article/ Numero di catalogo/Catalogue number
	Packungsinhalt/Contenu du coffret/ Contenuto della confezione/Contents of kit
	In vitro diagnosticum/Usage in vitro/ Per uso diagnostico in vitro/For in vitro diagnostic use
	Hersteller/Fabricant/Fabbricante/Manufacturer
	Lagertemperatur/Température de conservation/ Temperatura di conservazione/Storage temperature
	Verwendbar bis/Date de péremption/ Data di scadenza/Use before
	Chargenbezeichnung / Numéro du lot/ Codice del lotto/Lot number
	Gebrauchsanweisung beachten/Consulter le mode d'emploi/ Leggere le istruzioni per l'uso/Consult the instructions for use
	Einmalgebrauch/Usage unique/Per uso unico/single use
	CE-Kennzeichnung/CE marquage/ Marchio CE/CE marking

Field of application


Transport medium, bacterial count and selective propagation of bacteria in urine.

Package contents

- 10 UrinAX dip inoculum culture media
- 10 patient labels

Principle

UrinAX consists of a sealed container with 3 culture media:

The CLED culture medium is used to determine the total bacterial count in the urine. The MacConkey Agar is used for the selective propagation of Enterobacteriaceae. The Enterococcus Agar is used for selective propagation of Enterococcus.

Composition

Basic composition in g/l

CLED Agar	MacConkey Agar
Peptone 8.0	Peptone 20
Meat extract 3.0	Lactose 10
Lactose 10	Bile salts 1.5
L-cystein 0.128	Neutral red 0.03
Bromthymol blue 0.02	pH 7.0
pH 7.5	
Enterococcus Agar	
Peptone 8.0	Esculin 1.0
Bile salts 20	pH 7.1
Ammonium iron(III) citrate 0.5	

Storage

UrinAX should be stored at 15–25 °C. Temperature fluctuations and drafts should be avoided, and incubation should occur at 35–37 °C within 24 hours. Do not freeze. Do not use the product after the expiry date specified on the package.

Collection of urinary specimens

The urinary specimen must be inoculated into the UrinAX culture medium immediately after collection. Test results can be affected if the patient has been treated with antibiotics. In this case, the test must be performed no sooner than 48 hours after taking the last dose of medication.

Performing the test

1. Unscrew the UrinAX tube and remove the culture medium without touching the cultures.
2. Dip the culture medium into the urine mid-stream until the agar surfaces are completely covered. If the urine quantity is insufficient, pour the urine onto the agar surfaces. Allow surplus urine to run off the culture medium. Dab off the last drops of urine with filter paper.
3. Return the culture medium to the tube and close the lid.
4. Complete the patient labels and affix them to the tube.
5. Incubate the tube in an upright position in an incubator at 37 °C for 16–24 hours.

In the case of a negative result, the cultures can be incubated for a further 24 hours to detect slow-growing bacteria.


Transportation of the inoculated UrinAX

The inoculated culture medium can be transported at 15–25 °C. Incubation should occur within 24 hours.

Interpretation of results

After incubation, assess the bacteria growth on CLED Agar. Compare the bacterial count on the CLED Agar side with the following readout:

Bacterial count on CLED Agar/ml

Notes: Bacterial counts $< 10^4</math> bacteria/ml are usually not considered pathological. Bacterial counts between $10^4</math> and $10^5</math> are an uncertain result and should be repeated. Bacterial counts of $> 10^5</math> bacteria/ml indicate a urinary tract infection.$$$$

A mixture of different bacteria strains on the UrinAX is very probably due to contamination of the urinary specimen.

MacConkey Agar:

Assess the presence of colonies. With a positive lactose reaction, the colonies are red. With a negative lactose reaction, the colonies are colorless. Carry out biochemical or immunological identification of the bacteria.

Enterococcus Agar:

Assess the presence of colonies. The growth of enterococcus is indicated by a dark-brown/blackish precipitate caused by esculin hydrolysis.

Quality control

Quality controls are carried out during the manufacture of every UrinAX lot. Growth characteristics of the medium can be tested with the following strains:

- E. coli ATCC 25922
- Staphylococcus aureus ATCC 25923
- Proteus mirabilis ATCC 12453
- Enterococcus faecalis ATCC 29212

Interpret the results after a 16 to 24-hour incubation period as follows:

E. coli ATCC 25922
Yellow growth on CLED medium and pink colonies on MacConkey.

Staph. aureus ATCC 25923
Growth on the CLED medium only. Yellow colonies.

Proteus mirabilis ATCC 12453
Blue colouration on the CLED medium and colorless colonies on the MacConkey Agar.

Enterococcus faecalis ATCC 29212
The enterococci form yellowish colonies on the CLED medium and brown to black colonies on the E medium.

Disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

References

NCCLS Publication: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Volume 24:14, 2004.